


Product Specifications (1/2)

Product Name	Pearl Barley Syrup	
Seller	KYOTO GRAIN SYSTEM CO., LTD. Shiseido Kyoto Bldg. 2F, 480 Kitafudodo-cho Kizuyabashi-sagaru Aburanokoji-dori, Shimogyo-ku, Kyoto city, Kyoto 600-8233 Japan TEL : +81-75-353-8833 FAX : +81-75-353-8837	
Manufacturer	KYOTO GRAIN SYSTEM CO., LTD. Ishikawa factory 311-2 Koto Kahokugun Tsubatamachi, Ishikawa 929-0328 Japan TEL : +81-76-289-7794 FAX : +81-76-289-7796	
Material	Raw Material	Pearl Barley(dehulled), Malt
	Country Of Origin	Laos, Canada
	Allergens	No
	GMO	No
	BSE	None
Storage	Temperature	Room temperature
	Description	Store at cool and dark place avoiding hot and humid condition and exposure to strong light and direct sunlight.
Shelf life	Before opening	1080 days Production day (packing day) + 1080 days
Packaging	External	Cardboard box 300x300x300H mm 
	Internal	Plastic bag 300x300x300H mm 
Net Weight	28kg	
Code/Label		
Manufacturing process	<p>【Steaming】 Material input → soaking → Steaming → Temporary Storage</p> <p>【Sweetening】 Input → liquefaction/Saccharification → Filtration → Evapolation → Filtration → Filling → measuring → Packing → Labeling → Storage → Shipment</p>	

Product Specifications (2/2)

Product Name	Pearl Barley Syrup
---------------------	---------------------------

Product characteristics

Items	Standard value	Measuring method
Appearance	Clear, Gold Color	Visual inspection
Aroma & Taste	Toasty aroma without off-flavor and off-taste.	Sensory test
Foreign materials	Free from any unusual foreign materials.	Visual inspection
Brix	Over 80	Brix analyzer
		Particle size analyzer

Items	Unit	Standard value	Result	Measuring method
Arsenic as AS ₂ O ₃	ppm	-	-	-
Heavy Metal as Pb	ppm	-	-	-
general viable bacteria	cfu/g	Max 50,000	Meet with standard	Standard agar plate, incubation at 35°C for 48hours
Coliforms		Negative	Meet with standard	3M Petrifilm Coliform Count Plate, incubation at 35°C for 24hours
Mould, Yeast	cfu/g	Max 1,000	Meet with standard	Potato dextrose plate Count, incubation at 25°C for 7days
Staphylococcus aureus	cfu/g	Negative	Meet with standard	Salt Egg Yolk Agar Base
Pesticide residue		-	-	-
				-